

01-110X
Vortex Mixer

The Vortex Mixer is suitable for a variety of applications, from gentle mixing to vigorous resuspension of biological and chemical components in liquids. It is widely used in the fields of life science and physical and chemical analysis.



Features

- Orbital mixing movement
- Speed range: 0 - 2500 rpm, stepless speed settings
- Touch and continuous operation settings
- Cast aluminum base and specially designed vacuum suction feet keep the mixer stable
- Various adapters and accessories available for use with tubes of different sizes (sold separately)

