

01-220X Levo Plus Pipette Filler

The Levo Plus pipette filler efficiently aspirates and dispenses liquids, its powerful motor filling a 25 mL pipette in less than 5 seconds. It is compatible with most standard glass and plastic pipettes with volumes ranging from 0.1 - 100 mL.

Features

- Lightweight, ergonomic design
- Convenient one-handed operation
- Eight aspiration and dispensing speeds
- LCD display shows speed settings and remaining battery life
- Efficient lithium-ion battery can last for up to 8 hours of intermittent use
- Replaceable hydrophobic 0.45 µm filter
- Integral, autoclavable silicone adapter suits use with 0.1 - 100 mL pipettes
- Can be used while charging



Specifications

| | |
|--------------------------|--|
| Aspiration Speeds | 8 |
| Dispensing Speeds | 8 |
| Battery | lithium-ion |
| Battery Life | 8+ hours of intermittent use |
| Charging Time | 2-3 hours |
| Compatible Pipette Types | glass or plastic pipettes (0.1-100mL), Pasteur pipettes |
| Filter | 0.45 µm hydrophobic |