

## 01-221X Levo ME Pipette Filler

**NEW**

The Levo ME Pipette Filler has improved ergonomics and streamlined function for a more convenient pipetting experience.

### Features

- Ergonomic shape and intuitive operation
- Pressure-sensitive buttons: increasing pressure increases aspiration/dispensing speed
- Autoclavable hydrophobic filter prevents over-aspiration
- Rechargeable, long-lasting battery
- LED light indicates aspiration/dispensing, charging, and low battery
- Compatible with almost all standard glass or plastic pipettes



### Specifications

Catalogue No.	01-221X
Aspiration Speeds	aspirates 25mL of water in < 8 seconds
Dispensing Speeds	gravity dispensing
Battery	lithium-ion
Battery Life	8+ hours of intermittent use
Charging Time	2-3 hours
Compatible Pipette Types	glass or plastic pipettes (0.1-100mL), Pasteur pipettes
Filter	0.45 µm hydrophobic