

01-320X FlatSpin Magnetic Stirrer

This ultra-flat magnetic stirrer is easily portable and saves precious space, making it an excellent addition to any laboratory.

Features

- Ultra-flat, space-saving design
- Stirring speed range: 15 - 1500 rpm
- Automatic rotation-reversing stirring for better mixing results
- Constructed of chemically resistant materials
- Slip-resistant and stable



Specifications

Catalogue No.	01-320X
Voltage [VAC]	100-240
Frequency [Hz]	50/60
Motor input rating [W]	5
Motor output rating [W]	3
Stirring speed range [rpm]	15-1500
Max. stirring quantity [L, H ₂ O]	0.8
Max. load capacity [kg]	2.0
Dimension of work plate [mm]	90 × 90
Material	ABS
Dimension [W × D × H, mm]	175 × 120 × 15
Weight [kg]	0.3
Permissible ambient temperature [°C]	5-40
Permissible relative humidity	80%
Protection class according to DIN EN60529	IP65

