

03-4221, 03-4222 Elite Series Dry Bath Incubators

The Elite Series Dry Bath Incubators feature digital microprocessor control, user temperature calibration, and user-friendly design. One-block and dual-block types are available.

Features

- PID temperature control maintains temperature accuracy
- Leak-proof molded aluminum alloy chamber with PTFE coating
- Multi-functional: can also be used as a water bath or bead bath
- Various interchangeable heating blocks available




03-4221
One Block Unit



03-4222
Dual Block Unit

Dry Bath Blocks

These precisely machined aluminum alloy dry bath blocks deliver efficient heat transfer and are suitable for use with tubes ranging from 0.2 to 50 mL.



Specifications

Catalogue No.	03-4221	03-4222
Display	LCD Display	
Temperature Control Range [°C]	5 above ambient to 150	
Temperature Increment [°C]	0.10	
Temperature Accuracy at 37°C [°C]	± 0.20	
Operation modes	continuous or timed up to 99 hr:59 min	
Dry Bath Blocks	Capacity	
03-5002	0.2mL tubes x 64 or PCR 8-strip tubes x 8	
03-5005	0.5mL tubes x 20	
03-5015	1.5mL or 2.0mL tubes x 20	
03-5017	15mL centrifuge tubes x 12	
03-5029	50mL centrifuge tubes x 4	