

### 03-6200

### BluView Transilluminator

The BluView Transilluminator provides an outstanding gel viewing experience, utilizing a safe blue light that will not degrade or damage samples or pose a hazard to users.

- Sturdy aluminum alloy casing
- 16 x 16 cm gel viewing surface
- Lightweight, portable design



#### Specifications

| Catalogue No.                          | 03-6200            |
|--|--------------------|
| Blue Light Wavelength [nm]             | 470                |
| Illuminator Base Design                | Flat bed           |
| Construction Material                  | Aluminum alloy     |
| Dimensions [W x L x H, mm]             | 200 x 200 x 28.7   |
| Weight [kg]                            | ~1.5               |
| Power Adaptor                          | 12VDC/5A, 60W Max. |
| Filter Viewing Surface Dimensions [cm] | 16 x 16            |