

01-120X, 01-130X Classic Rotators

These classic adjustable-speed rotators provide continuous end-over-end rolling action to mix samples thoroughly or keep them in suspension. Two types are available.

Features

- Speed range: 0 - 80 rpm
- End-over-end rotation
- Continuous operation
- Compact, space saving design
- Various adapters available for use with 1.5mL - 50mL tubes



01-120X (radial)



01-130X (axial)

Specifications

Catalogue No.	01-120X, 01-130X
Voltage [VAC]	100-240
Frequency [Hz]	50/60
Power [W]	20
Motor type	DC motor
Tilt angle	90°
Speed range [rpm]	0-80
Operation type	Continuous
Dimensions [D×W×H, mm]	300 x 220 x 310 (01-120X) 150 x 530 x 190 (01-130X)
Weight [kg]	2.0 (01-120X) 1.7 (01-130X)
Permissible ambient temperature [°C]	5-40
Permissible relative humidity	80%
Protection class according to DIN EN60529	IP21

