

01-800X Palm Micro Centrifuge

The Palm Micro Centrifuge is a very small, compact centrifuge ideal for microfiltration and quick spins. It has been tested and approved in accordance with IEC/EN61010-2-20.

Features

- Fixed speed: 7000 rpm
- Quick acceleration and deceleration
- Rotors for microcentrifuge tubes, PCR tubes, and PCR 8-strip tubes included
- Easily change rotors without tools
- Quiet, low-noise design



CE FC



Specifications

Catalogue No.	01-800X
Max. speed [rpm]	7000
Max. RCF	2680 x g
Rotor Capacity	0.2/0.5/1.5/2.0mL x 8, or 0.2mL PCR 8-strips x 2
Run Time	Continuous operation
Driving Motor	DC motor
Power	AC 110-240V, 50Hz/60Hz
Noise level [dB]	≤45
Dimensions [D x W x H, mm]	150 x 150 x 117
Weight [kg]	0.5